

## Upholstery Colours

The following colours of high quality synthetic leather are available for the upholstery of our products. The synthetic leathers we use are characterized by their tensile strength, durability, longevity and easy cleaning and resistance to disinfectants. We use a colour selection of synthetic leathers VOWAled Pisa (Vowalon) and MEDISTAT (Cotting/Griffine).



The synthetic leather should be cleaned with a soft brush and a mild soap solution (diluted commercial dishwashing liquid - no organic solvents). After cleaning, always rinse with clear water and dry with a soft cloth.

We also offer the option of having the upholstery done in other colours or brands of synthetic leather. The corresponding materials should be sent to us, subject to prior arrangement.

The black electrically conductive synthetic leather upholstery (with earth contact on the frame) is available at extra charge. It is, however, already included in the price of the operation and multi-function tables.