

## Battery powered X-Ray tables AGA-AKKU-POWER-LIFT and AGA-AKKU-POWER-DUO-LIFT

The X-ray tables AGA-POWER-LIFT (with electric height adjustment) and AGA-POWER-DUO-LIFT (with electric height and tilt adjustments) have been specially developed for use in the field of radiology. The absence of a longitudinal brace at the front of the leg frame makes the use of a C-arm easily possible.

These two X-ray tables are now available in a battery powered version. The new X-ray table models presented below, AGA-AKKU-POWER-LIFT (with height adjustment) and AGA-AKKU-POWER-DUO-LIFT (height and tilt adjustments) are equipped with rechargeable batteries for greater versatility of use.

The complete product descriptions with additional features, as well as the available options and accessories may be found in our new catalogue (issue 42) on pages 24-27.



1 RGL-1065/NAC, PAL-190-65, LR/Z-125 (H: 620-1020 mm)



2 RGL-1065/NAC, LR/Z-125 (H: 620-1020 mm)

	AGA-AKKU-POWER-LIFT		AGA-AKKU-POWER-DUO-LIFT	
Article number	RGL-1065/AC	RGL-1080/AC	RGL-1065/NAC	RGL-1080/NAC
Picture(s)	Not shown	Not shown	In 1, 2	Similar 1, 2
Price (EUR)				
Length (mm)	1900		1900	
Tabletop	Continuous		Continuous	
Tabletop in 2 sections (500/1400 mm)	RGL-1065/KAC	RGL-1080/KAC	RGL-1065/NKAC	RGL-1080/NKAC
Width (mm)	650	800	650	800
Height adjustment (mm)	460-860		460-860	
Tilt adjustment	_		17°	
Height adjustment mechanism	Battery powered, with hand switch, 2 buttons (up/down)		Battery powered, with hand switch, 6 buttons for simultaneous height and tilt adjustment	
Maximum load (kg)	150		150	
Battery module	BMEAC		BMEAC	
Paper roll holder	P-2065/A	P-2080/A	P-2065/A	P-2080/A
Accessories	Catalogue page 27		Catalogue page 27	