

(OR-) Swivel Stool with foot release

We have extended the range of our swivel stools (with hand lever release) to include models with foot release. The swivel stools with foot release, when equipped with electrically conductive double castors and synthetic leather cover, are particularly well suited for use in the operating room.



1 DHFA-2000/LEIT



2 DHFA-2000

OR Swivel Stool DHFA-2000/LEIT

- Five-arm aluminium foot base (Ø 610 mm) in a polished finish
- Electrically conductive braked soft double castors
- Seat Ø 380 mm, upholstery thickness 75 mm
- Seat upholstered with electrically conductive synthetic leather (colour: black).
- Height adjustment of the seat (approx. 460-580 mm) activated with the foot via a foot ring located under the base
- Chrome-plated foot ring (Ø 300 mm)
- Maximum load 140 kg
- Product weight approx. 9 kg

Swivel Stool DHFA-2000

- Five-arm aluminium foot base (Ø 610 mm) in a polished finish
- Braked (hard) double castors for soft floors.
- Seat Ø 380 mm, upholstery thickness 75 mm
- Seat upholstered with synthetic leather (21 colours available)
- Height adjustment of the seat (approx. 460-580 mm) activated with the foot via a located under the base.
- Chrome-plated foot ring (Ø 300 mm)
- Maximum load 140 kg
- · Product weight approx. 9 kg

| Picture | Article Number | Adjustment Range (mm) | Aluminium Base Polished | Price (EUR) |
|---------|----------------|--------------------------|-------------------------|----------------|
| 1 | DHFA-2000/LEIT | 460-580 | х | |
| 2 | DHFA-2000 | 460-580 | x | |