

(OR) Swivel Chairs with foot release

We have extended the range of our swivel chairs (with hand lever release) to include models with foot release. The swivel chairs with foot release, when equipped with electrically conductive double castors and synthetic leather cover, are particularly well suited for use in the operating room.



1 DSFA-2000/RLLEIT



2 DSFA-2000/RL

OR Swivel Chair DSFA-2000/RLLEIT

- Five-arm aluminium foot base (Ø 610 mm) in a polished finish
- Electrically conductive braked soft double castors
- Seat Ø 380 mm, upholstery thickness 75 mm
- Seat and backrest upholstered with electrically conductive synthetic leather (colour: black).
- Slight inclination (3°) of seat towards the backrest
- Adjustable ergonomic backrest
- Height adjustment of the seat (approx. 490-610 mm) activated with the foot via a chrome-plated foot ring (Ø 300 mm) located under the base.
- Maximum load 140 kg, product weight approx. 10 kg

Swivel Chair DSFA-2000/RL

- Five-arm aluminium foot base (Ø 610 mm) in a polished finish
- Braked (hard) double castors for soft floors.
- Seat Ø 380 mm, upholstery thickness 75 mm
- Seat and backrest upholstered with synthetic leather (21 colours available)
- Slight inclination (3°) of seat towards the backrest
- Adjustable ergonomic backrest
- Height adjustment of the seat (approx. 490-610 mm) activated with the foot via a chrome-plated foot ring (Ø 300 mm) located under the base.
- Maximum load 140 kg, product weight approx. 10 kg

Picture	Article Number	Adjustment Range (mm)	Aluminium Base		Price (EUR)
			Polished		
1	DSFA-2000/RLLEIT	490-610	x		
2	DSFA-2000/RL	490-610	x		