

## Privacy Screens and Miscellaneous Furniture

---







# Privacy Screens

## Privacy Screens

The folding privacy screens, also referred to as room dividers, can be considered indispensable pieces of equipment in every medical practice, clinic or hospital.

They are available in three versions: stationary without castors, mobile with castors or pull-out/wall-mounted with castors. The sturdy frame is made of profiled anodized aluminium tubing (30x20 mm) which offers the advantage of a very low weight. The individual panels are joined by two galvanized hinges that can be folded in both directions.

For safety and best stability, the panels should be set up at an angle of 135° at all times.



1 PV-3004



2 PV-3005

Privacy Screens				
Article number	PV-3002	PV-3003	PV-3004	PV-3005
Picture(s)	3	4	1	2
Price				
Height (mm)	1652	1652	1652	1652
Number of panels	2	3	4	5
Weight (kg)	7,70	11,60	15,50	19,40
Dimensions (mm)				
▪ Extended (length)	1130	1690	2250	2810
▪ Folded (LxD)	605x56	605x82	605x105	605x120



3 PV-3002



4 PV-3003



5 Colour options for privacy screen panels

# Privacy Screens

The panels made of 3 mm thick, impact resistant and easy-to-clean PVC rigid foam boards. The even, micro-cell foam structure is uniformly coloured, flame retardant and resistant to repeated disinfectant use.

The privacy screen panels are available in several colours which may be combined as desired, at no extra charge. Please specify your selection upon ordering.



6 PV-3005/LR



7 PV-3003/LRWM

	Mobile Privacy Screens				Wall-Mounted Privacy Screens			
Article number	PV-3002/LR	PV-3003/LR	PV-3004/LR	PV-3005/LR	PV-3002/LRWM	PV-3003/LRWM	PV-3004/LRWM	PV-3005/LRWM
Picture(s)	Not shown	Not shown	8	6	Not shown	7	Not shown	Not shown
Price								
Height with castors (mm)	1745	1745	1745	1745	1745	1745	1745	1745
Number of panels	2	3	4	5	2,5	3,5	4,5	5,5
Weight (kg)	8,65	12,65	16,65	20,65	11,40	15,40	19,40	23,40
Dimensions (mm) • Extended (length) • Folded (LxD)	1130 605x56	1690 605x82	2250 605x105	2810 605x120	1410 605x92	1970 605x118	2530 605x131	3090 605x146
Castors (Ø 65 mm)	✓	✓	✓	✓	✓	✓	✓	✓

## Mobile Privacy Screens

- The bottom extremities of the two outside panels are each fitted with a rotating aluminium crossbar. Each crossbar has one locking and one non-locking double castor.
- The inner panels are also fitted at their junction with one non-locking double castor.

## Wall-Mounted Privacy Screens on Castors

- An aluminium rail with screws and dowels for wall-mounting are included with each product.
- The bottom of the outside panel is fitted with a rotating aluminium crossbar.



8 PV-3004/LR

# Miscellaneous Furniture

## Step Stools

The step stools are available in two frame options (chrome-plated or stainless steel) and with one or two steps. The frame is made of square tubing profiles (20x20 mm). The slip-resistant steps are made of black polyurethane (PU) coated foam plates.



1 FA-3000/1



2 FA-2000/2


Article Number	FA-2000/1	FA-2000/2	FA-3000/1	FA-3000/2
Picture(s)	similar 1	2	1	similar 2
Price				
One step	✓	—	✓	—
Two steps	—	✓	—	✓
Step height(s) (mm)	240	220 440	240	220 440
Frame	Chrome		Stainless Steel	
Maximum load (kg)	150		150	
Dimensions of the plates (mm)	490x320x15 (WxDxH)		490x320x15 (WxDxH)	

## Infusion (IV) Holders

The infusion holders are available in two versions: mobile on castors and wall-mounted. The mobile infusion holder may be ordered under Article No. INF-20215. The infusion holder (Article No. IWH-2045) is intended for wall-mounting.



3 INF-20215



4 IWH-2045

Article Number	INF-20215	IWH-2045
Picture(s)	3	4
Price		
Length (mm)	1350-2150	505 (620 with cup)
Five-arm base with double castors (of which two are conductive and locking)	Aluminium Ø 630 mm	—
Frame	Chrome	Chrome
Bottle hooks	4	2
Max. load each hook (kg)	2	2
Height adjustment range of the hooks (mm), secured by screw lock mechanism	1350-2150	—
Plastic drip cup and holder	✓	✓
Castors (Ø 50 mm)	✓	—
Fastening material for wall-mounting	—	Dowels and screws

## Miscellaneous Furniture

### Drawer Units and Multi-Purpose Mobile Cart

The drawer units, mobile drawer units and multi-purpose mobile cart offer versatile, mobile storage for supplies and equipment. The body of the units and the drawer fronts are made of white melamine-coated panels, thickness 16 mm.



1 SW-470/2, SW-470/3, U-1065



2 UF-580/2



3 HB-580/3MBG

Article Number	SW-470/2	SW-470/3	UF-580/2	UF-580/3	HB-580/2	HB-580/3
Picture	In 1	In 1	2	Similar 2	Similar 3	Similar 3
Price						
Dimensions (mm) incl. fronts, handle and castors	600x540x550 (WxDxH)		400x538x570 (WxDxH)		400x538x570 (WxDxH)	
Width of drawer fronts (mm)	557		433		433	
Castors (Ø 50 mm, white)	✓		✓		✓	
Two drawers	✓	–	✓	–	✓	–
Three drawers	–	✓	–	✓	–	✓
Padded seat and front upholstered with the same synthetic leather as the table	–		✓		–	
Drawer fronts with bare handles	✓		–		HB-580/2MBG	HB-580/3MBG
Raised edge (10 mm) on three sides	✓		–		✓	

### Multi-Purpose Mobile Cart

Article Number	VZW-1062
Picture	4
Price	
Frame	Lacquered
Castors (Ø 75 mm)	✓
Dimensions (mm) incl. 30 mm deep handle	620x530x735 (WxDxH)
Dimensions of shelves (mm)	580x460 (WxD)
Height of shelves (mm)	225/665
Drawer with chrome-plated handle	✓

#### Additional Features

- The frame of the multi-purpose mobile cart is made of welded rectangular steel pipe, lacquered in pure white RAL9010 (unless another colour is specified at time of order, see p. 4).
- The multi-purpose mobile cart is delivered fully assembled, with only the four black nylon double castors to be inserted in the corresponding sleeves.



4 VZW-1062