

# AGA

## Sanitätsartikel

### Declaration of Compatibility for MRI products

Manufacturer: AGA Sanitätsartikel GmbH  
Löhner Strasse 198b  
D-32584 Löhne, Germany



Products: MRI Patient transport tables

We hereby declare that the above mentioned medical products were independently tested and are compatible with MRI systems **up to 7 Tesla**.

Article Numbers: U-3065/LRVA, EKG-3080/LRVA, AWL-3065/VA, AWL-3080/VA

The MRI patient transport tables can be used safely under the following conditions :

- static magnetic field of maximum 7 Tesla or less,
- with an additional spatial gradient of maximum 7 T/m or less,
- the table must not exceed the limit line of 7 T/m with the head end/foot end in the direction of movement towards the magnet.

Article Numbers: H-U 3065/Y 980LRVA, H-EKG 3080/Y 980LRVA, H-AWL 3065/Y 980VA, H-AWL 3080/Y 980VA

The MRI patient transport tables can be used safely under the following conditions :

- static magnetic field of maximum 7 Tesla or less,
- with an additional spatial gradient of maximum 10 T/m or less,
- the table must not exceed the limit line of 10 T/m with the head end/foot end in the direction of movement towards the magnet.

Date: November 8th, 2023

Signature: Achim Walter  
Managing Director

